




Transforming Lives Through
Comprehensive Care




Ease of Use of AlignRT in DIBH Treatment

Dr Prithvi Raj Jampana MRCP FRCR
Consultant and Head, Radiation Oncology
Sindhu Hospitals, Hyderabad


Importance of DIBH in Breast Cancer

- Reduces radiotherapy dose to the heart and Lungs
 - Darby et al:
 - 2168 women had RT for breast cancer, 1205 controls, 963 events
 - Mean dose to heart: 6.6Gy for left breast and 2.9Gy for right breast, 4.9Gy
 - Rate of coronary event increased by 7.4% for each 1Gy increase mean dose
- | | |
|---------------------------|--|
| -% increase for mean dose | -% increase no of yrs since radiotherapy |
| 2Gy-----10% | 0-4 yrs-----16.3% |
| 2-4Gy-----30% | 5-9 yrs-----15.5% |
| 5-9Gy-----40% | 10-19yrs-----1.2% |
| 10Gy or more-----116% | 20 or more-----8.2% |



Age at irradiation (years)	Mean heart dose (Gy)	Cumulative risk (%) by attained age				Absolute risk (%) of radiation-related RHD death by age 80 years	Cumulative risk (%) by attained age				Absolute risk (%) of radiation-related RHD death by age 80 years
		50	60	70	80		50	60	70	80	
40	0	0.03	0.1	0.5	2.0	0.0	0.08	0.3	0.9	3.3	0.0
	0.5	0.03	0.2	0.5	2.0	0.0	0.09	0.3	0.9	3.4	0.1
	1	0.03	0.2	0.6	2.1	0.1	0.09	0.3	1.0	3.5	0.2
	2	0.03	0.2	0.6	2.2	0.2	0.10	0.3	1.0	3.8	0.5
	3	0.03	0.2	0.6	2.4	0.4	0.10	0.3	1.1	4.0	0.7
	4	0.03	0.2	0.7	2.5	0.5	0.11	0.4	1.2	4.2	0.9
	5	0.04	0.2	0.7	2.7	0.7	0.11	0.4	1.2	4.5	1.2
	6	0.04	0.2	0.8	2.8	0.8	0.12	0.4	1.3	4.7	1.4
	7	0.04	0.2	0.8	3.0	1.0	0.13	0.4	1.4	4.9	1.6
	8	0.04	0.2	0.8	3.1	1.1	0.13	0.5	1.4	5.2	1.9
	9	0.04	0.2	0.9	3.2	1.2	0.14	0.5	1.5	5.4	2.1
50	10	0.05	0.3	0.9	3.4	1.4	0.15	0.5	1.6	5.6	2.3
	0	0.09	0.5	1.9	0.0	0.3	0.9	3.4	0.0		0.0
	0.5	0.09	0.5	2.0	0.1	0.3	1.0	3.5	0.1		0.1
	1	0.09	0.5	2.1	0.2	0.3	1.0	3.6	0.2		0.2
	2	0.10	0.5	2.2	0.3	0.3	1.1	3.8	0.4		0.4
	3	0.11	0.6	2.4	0.5	0.3	1.1	4.1	0.7		0.7
	4	0.11	0.6	2.5	0.6	0.4	1.2	4.3	0.9		0.9
	5	0.12	0.6	2.7	0.8	0.4	1.3	4.6	1.2		1.2
	6	0.12	0.7	2.8	0.9	0.4	1.3	4.8	1.4		1.4
	7	0.13	0.7	2.9	1.0	0.4	1.4	5.0	1.6		1.6
	8	0.14	0.7	3.1	1.2	0.4	1.5	5.3	1.9		1.9



Typical DIBH options in Breast Cancer




- Voluntary breath hold
 - UK Heartspare Study—Radiother Oncol. 2013 Aug;108(2):242-7. doi: 10.1016/j.radonc.2013.04.021. Epub 2013 May 29.
-<https://pmc.ncbi.nlm.nih.gov/articles/PMC4211647/>-
 - Easy to use, no need for any equipment
 - More suitable with 3D-CRT, close monitoring during treatment
 - Outcomes similar but Patients felt more comfortable with vDIBH
- Spirometry based Breath hold
 - Tracks patient's lung volume, needs cooperation from patient
- Surface based systems
 - Use projectors and cameras for real time 3d surface of the chest
Eg: SGRT



Our DIBH workflow using SGRT


- Patient explained to about DIBH- suitability assessed
- Patient seen in RT department-explained about DIBH procedure
- Ideally, at least 1 day time given to practice breath hold
 - up to 80% of full inspiratory capacity
 - check chest vs abdominal breathing
- Patient taken for planning CT scan
- Mirror placed on CT scan allows CT technicians to observe breath hold consistency



Our DIBH workflow using SGRT



- Patient practices breath hold on CT Couch
- Skin marks placed using lasers in free breathing and breath hold positions
- Displacement of skin marks noted between free breathing and breath hold positions
>6mm, we consider suitable, 1cm-1.2cm, ideal
- CT scan done in free breathing and Breath hold



Our DIBH workflow using SGRT



- Treatment planned with 3D-CRT using field in field technique, IMRT or VMAT
- If VMAT, we plan with 2 tangential arcs except when IMN also treated

• Heart doses aimed for while planning	40Gy/15#	26Gy/5#
	V10 < 5%	V7Gy < 5%
	V2 < 30%	V1.5Gy < 30%

- Mean Heart Dose-----3.2Gy, max 4.0Gy (RTOG 1005)
- Our Data for Mean heart dose----2.29Gy with DIBH and 3.6Gy without DIBH

Our DIBH workflow using SGRT

- Patient set up in position using postural video, position verified in free breathing and breath hold using SGRT
- CBCT done in breath hold to confirm accurate positioning
- Once CT verified, patient instructed to breathe in and hold breath, beam switches on. Breath hold held for 20 secs.
- Patient specifically instructed to breathe in and hold for first 2-3 fractions
- After that, most patients able to do the procedure by themselves




Our DIBH workflow using SGRT



Our DIBH workflow using SGRT

- Chest Wall treatment with bolus



In Conclusion



- DIBH an essential part of treatment in breast cancers, especially left breast cancers
- Helps reduce heart dose and also lung doses
- Different techniques available to do DIBH
- SGRT using AlignRT, in our experience
 - Excellent patient comfort
 - Very accurate
 - Easy to setup and implement
 - Easy for RT Technologists to monitor patients
 - More consistent and measurable results

Thank You



Sindhu Health Campus, HITEC City
Hyderabad-500084, Telangana, India.
Ph: (040) 4545 4545
email: info@sindhuhospitals.com