

Simulation Free Radiotherapy (SimLess) Using Surface Guided Radiotherapy (SGRT)

Sally Evill

Radiation Therapist
GenesisCare Crowsnest Mater


SGRT Annual APAC Meeting 2026





What is SimLess?

SimLess RT allows a patient to be planned and treated without the need for a specific CT appointment performed in the RT department.



The Challenge

Setup Reproducibility in SimLess RT

- Diagnostic imaging \neq Treatment position
- Limited immobilisation
- Increased uncertainty in setup
- Potential need for reimaging





SGRT in Our Department

01 AlignRT used as a standard of Care

02 Routine for setup and monitoring

03 High Staff familiarity

04 Improves positioning confidence

05 Intrafractional monitoring

06 Reduces set up time

07 Potential in decreasing re-imaging

08 Real-time surface Matching

AlignRT Preparation for SimLess Treatment

Import Structures

- Approved Plans Structure set sent to AlignRT

Select Protocol



- Choose Protocol (Usually Pelvis)
- Tolerance 0.5cm

Define ROI

- Region of interest set for patient set up


Verification

- Check during final plan & Pre-treatment check



AlignRT Workflow on Treatment

- ⑦ AlignRT integrated with Linac
- ⑦ Set-up with video and Postural Alignment features
- ⑦ Full CBCT taken for every treatment
- ⑦ Re-image is shifts larger then 2cm
- ⑦ Recapture taken and Monitoring enabled




The software interface displays patient alignment data and a 3D model of the pelvis. The patient is identified as "SGRT Patient" and the treatment area is "Pelvis". The alignment data is as follows:

| Parameter | Value |
|--------------------|-------|
| VRT _{cm} | -0.12 |
| LNG _{cm} | -0.02 |
| LAT _{cm} | 0.11 |
| MAG _{cm} | 0.17 |
| RTN [°] | -0.3 |
| ROLL [°] | 0.9 |
| PITCH [°] | -0.2 |

The 3D model shows the pelvis with a cyan surface and magenta contours. The couch angle is 0.0°. The beam control is ON. The interface also shows a "BEAM: ENABLED" warning and a "Coaching" dropdown menu with options for "Surface", "Deformation", and "Video".

Evaluation of Workflow

- ⑦ 10 palliative patients
- ⑦ Pelvic or Spine Bone Metastases
- ⑦ Planned using a recent Diagnostic Image
- ⑦ No dedicated CT simulation
- ⑦ All patients in room set up with AlignRT & Intrafractional Monitoring




Results- Feasibility & Efficiency

- ✓ 10/10 patients treated successfully using SimLess + SGRT
- ✓ Average treatment time 19 minutes (Within standard 20 min slot)
- ✓ No cases reverted to traditional CT simulation workflow (More than 10 CT hours Saved)
- ✓ Reduced need for reimaging
- ✓ Less than 3% PTV Dmax Difference

Results- Setup Accuracy

Translational Shifts



Retrospective Dose Analysis

- ✓ Physics team analysed dose evaluations on SimLess treated patients
- ✓ Original plans are calculated onto the treated CBCT using fixed MUs
- ✓ This showed a mean PTV difference of Less than 3%
- ✓ Further validating the SimLess technique

Clinical Impact



Positioning Confidence

Improved confidence in setup without dedicated planning CT



Efficiency

Delivered within standard 20-minute appointment slots



Streamlined Pathway

Simplified palliative workflow from planning to treatment



Patient & Staff Experience

Positive informal feedback from both groups

Limitations

- ① Small Cohort (n=10 patients)
- ① Single Centre Clinical Experience
- ① Limited Tumour Sites
- ① No direct comparison with CT-simulation Workflow

Supporting SimLess Radiotherapy Using SGRT

- ➔ SGRT integrates seamlessly into the SimLess pathway
- ➔ Accurate positioning achieved without planning CT
- ➔ Delivered within standard clinical timeframes
- ➔ Supports efficient, patient-centred palliative care



Questions and answers



Sally Evill

Speaker

Contact Details



Sally.evill@genesiscare.com



Phone number

Thank you

Lorem ipsum dolor sit amet, consectetur adipiscing elit, sed do eiusmod tempor incididunt ut labore et dolore magna aliqua. Duis tristique sollicitudin nibh sit amet commodo nulla. Elit eget gravida cum sociis. Varius morbi enim nunc faucibus a pellentesque sit amet.

